
MEDICA DEPOT

# BELOTERO® BALANCE

## Product Specifications

|    |                                            |
|----|--------------------------------------------|
| F  | NOTICE BELOTERO® BALANCE                   |
| GB | INSTRUCTIONS FOR USE FOR BELOTERO® BALANCE |
| D  | GEBRAUCHSINFORMATION FÜR BELOTERO® BALANCE |
| I  | ISTRUZIONI PER L'USO DI BELOTERO® BALANCE  |

Aillette en position adéquate pour l'injection  
Backstop in the right position during injection  
Korrekte Position der Fingerauflage bei der Injektion.  
Posizione corretta del pistone durante l'iniezione

|   |                          |
|---|--------------------------|
| F | NOTICE BELOTERO® BALANCE |
|---|--------------------------|

#### Description

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## GB INSTRUCTIONS FOR USE FOR BELOTERO® BALANCE

#### Description

BELOTERO Balance is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer.

#### Presentation

BELOTERO Balance is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

#### Composition

Cross-linked sodium hyaluronate: 22.5 mg/ml

Phosphate buffer pH 7 q.s.:

gel volume

The volume of the gel in each syringe is stated on the external box.

#### Indications

BELOTERO Balance is an injectable resorbable implant indicated to fill moderate facial wrinkles and folds as well as for lip enhancement. It is also suitable for correction of facial atrophic scars.

#### Posology and administration method

BELOTERO Balance is designed to be injected into the superficial to mid-dermis by a legally approved practitioner.

BELOTERO Balance can be used for all skin types.

For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for filling wrinkles and folds. A good knowledge of the anatomy and physiology of the site to be treated is required.

The treatment must be carried out under appropriate aseptic conditions.

BELOTERO Balance must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.

If necessary a local or loco-regional anaesthetic may be applied according to its instructions for use.

It is recommended to use one of the supplied needles.

To ensure optimal use of BELOTERO Balance, it is recommended to assemble the needle according to the diagrams below.

The practitioner should note that the use of 30G½ needles requires greater force to inject the gel than 27G½ needle.

BELOTERO Balance can be used in combination with other BELOTERO products during the same session. Instructions for use of each product should be followed.

If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.

BELOTERO Balance should be injected slowly. The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.

Gently massage the treated area after the injection to distribute the product uniformly.

#### Contra-indications

BELOTERO Balance is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate,
- In pregnant and breast-feeding women,
- In young patients under 18 years old,
- In patients presenting a general infection.

Do not inject BELOTERO Balance into blood vessels.

Do not inject BELOTERO Balance into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...).

Do not inject BELOTERO Balance into an area previously treated with a permanent dermal filler.

#### Precautions for use

Before treatment the patient must be informed about the device, its contra-indications and possible side effects.

In the absence of available clinical data on tolerance and efficacy of the injection of BELOTERO Balance in patients with antecedents or with an active autoimmune disease, or in patients presenting a history of severe mul-

tiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Balance on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject BELOTERO Balance in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or cheloids.

BELOTERO Balance can be injected in the glabellar area only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area could potentially lead to local vascular occlusion, ischemia and necrosis.

BELOTERO Balance can be injected in the periorbital region only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may be associated with an increased frequency and severity of side effects.

BELOTERO Balance must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO Balance must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO Balance with the above-mentioned treatments.

No clinical data is available on the injection of BELOTERO Balance into an area already treated with another filling product.

Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.

Patients using antithrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas, nodules or bleeding at the injection site.

Do not transfer BELOTERO Balance into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe.

BELOTERO Balance must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.

Discard the syringe and the remaining product after use.

Do not re-sterilize and do not reuse due to the associated risks including infection.

The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

#### Incompatibilities

Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Balance does not come into contact with such substances.

There is no known interaction with other local or loco-regional anesthetics.

#### Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema or pain sometimes accompanied by itching in the treated area. These reactions may last for a week.
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially to sodium hyaluronate.
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

#### Assembly of needle to syringe

For optimal use of BELOTERO Balance, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. **Firmly hold** the glass cylinder of the **syringe** and the **Luer-lock adaptor** between the thumb and the index fingers.
2. Grasp the protective cap with the other hand and unscrew it.
3. **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

#### Storage

Store between 2 °C and 30 °C. Protect from light and freezing. Avoid mechanical shocks.

#### References

Updated documentation may be available from ANTEIS SA, Switzerland.

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# MEDICA DEPOT

Medica Depot does not claim ownership of the content presented in this product insert - it is a photo reproduction of the product leaflet provided by the

Medica Depot team for your convenience. The intellectual property rights, including but not limited to copyrights, trademarks, and other applicable rights, associated with the content of this product insert remain the sole property of the original manufacturer.

Medica Depot disclaims any responsibility for any inaccuracies, omissions, or discrepancies in this document. Your use of this product insert is subject to our terms and conditions, available on our website. By accessing and using this document, you acknowledge and agree to the mentioned terms and conditions.