MEDICA DEPOT

HYACORP FINE®

Product Specifications

Skin Rejuvenation.

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INSTRUCTIONS FOR USE

Description

HYAcorp Fine is an absorbable skin implant with a high level of purity. It is a medical device produced from a hyaluronic acid that is not derived from animals. **HYAcorp Fine** is a sterile, apyrogenic, visco-elastic, biologically compatible (non-immunising, non-inflammatory, non-toxic) gel implant that is soluble in water and produced from a hyaluronic acid gained through fermentation. Hyaluronic acid is a naturally occurring polysaccharide in the dermal matrix of human skin. Hyaluronic acid is chemically, physically and biologically identical in the tissues of all higher organisms.

Composition

1 ml **HYAcorp Fine** contains:

Hyaluronic acid sodium salt	14.0 mg
Sodium chlorid	6.9 mg
Water for injection ad	1.0 ml

Mode of action

HYAcorp Fine is implanted into the medium dermal tissue to supplement the intercellular matrix and the intradermal tissue and to restore lost anatomical structures of the skin. Its mechanism of action is based on the latest biotechnology in the production of injectable hyaluronic acid.

Indication and application

- Replaces lost hyaluronic acid in the skin
- Increases elasticity
- Improves skin hydration
- Provides skin with a fresh look through its lifting effect

The result that can be achieved is dependent on skin type and the changes requested. The use is only permitted by medical personal.

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Contraindications

HYAcorp Fine may not be used in the event of:

- a tendency to hypertrophic and keloid scarring
- an intolerance towards gram-positive bacteria
- active inflammatory or infectious processes

- acute or chronic skin diseases
- anti-coagulant therapy
- a known allergy against hyaluronic acid

Patients with multiple allergies should be excluded from treatment. The use of **HYAcorp Fine** for breast augmentation is contraindicated.

NOTICE

HYAcorp Fine is intended for intradermal injection only and may not be injected into blood vessels as this could result in an occlusion of the vessels and an embolism. No clinical data are available on the application of **HYAcorp Fine** during pregnancy or lactation or to adolescents under 18 years of age. **HYAcorp Fine** is packaged sterile and is intended only for single use. It may not be re-sterilised. If the packaging is open/defective, it must not be used.

HYAcorp Fine is supplied in a sterile syringe and is ready-to-use. It may not be mixed with other injection agents.

Precautions

A complete anamnesis must be established prior to treatment to exclude possible contraindications. The normal precautions that apply to all intradermal injections must be observed. The implantation of **HYAcorp Fine** involves the risk of infection. The skin must be disinfected prior to treatment. Like other implants, **HYAcorp Fine** may not be applied near the site of inflammation or infection when treating patients with existing infectious or inflammatory processes. Before treatment, patients should not ingest aspirin, steroids or high doses of vitamin E as these substances may cause bleeding or raise susceptibility towards infection at the injection site. The area treated may not be exposed to great heat (sun, solarium, laser and IPL). The syringes and needles used are considered to be contaminated and must therefore be destroyed in accordance with the accepted rules of medical practice.

Adverse effects

Caused by the injection:

As with any other injection, patients may suffer from the following symptoms:

- temporary erythema
- slight swelling
- pain
- itching
- discolouration
- hardening

Typically these reactions spontaneously disappear within 2 to 5 days following the injection.

Caused by the product:

A hypersensitivity to hyaluronic acid following the injection has been reported less than 1 % out of 5000 treatments. This hypersensitivity is shown through an extended erythema, swelling and hardening at the implantation site. These reactions can occur immediately after the injection or up to 2 to 4 weeks later.

Application instructions

The areas to be treated must be marked prior to treatment. The treating person should note the facial expression and existing lines and a possible facial asymmetry. To perform the implantation as painfree as possible, a local anaesthetic should be used. The treating person must advise the patient of all precautions and possible adverse effects prior to treatment. The area to be treated must carefully be cleaned with antiseptic agents prior to treatment. The syringe is removed from the blister pack, the cap covering the tip of the syringe is removed and a needle is applied to the syringe. **HYAcorp Fine** is injected with the aid of 30G nanoneedles. The implantation is effected in the dermis.

Implantation technique

Implant **HYAcorp Fine** flatly in the medium dermis with the aid of the 30G nanoneedle. After implanting **HYAcorp Fine**, massage the treated area with slight pressure to guarantee the even distribution of the product. **HYAcorp Fine** sports excellent flowability in tissue.

ATTENTION

The graduation on the syringe is intended as an orientation aid for users with regard to final volume. It does not have any measuring function, it solely provides the amount used in relation to the nominal volume of 1 ml. Injection of a sufficient amount of the material is checked by the treating person visually and through touch.

Nature and contents of container

HYAcorp Fine is supplied sterile in a 1 ml syringe with an integrated Luerlock adapter in a blister pack for single use only. Sterile 30G nanoneedles, instructions for use and labels containing the lot number and the shelf life are packaged in a carton together with the blister pack. The patient is given one of the labels to guarantee traceability of the product.

NOTE

If the packaging and the blister pack are damaged the product must not be used. Do not re-sterilise **HYAcorp Fine**.

Precautions for storage

Do not freeze or expose to extreme heat.

HYAcorp Fine should be stored at room temperature (2 °C – 25 °C).

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CE 2409
BS_HYAcorp_Fine_BPZ_004

MEDICA DEPOT

Medica Depot does not claim ownership of the content presented in this product insert - it is a photo reproduction of the product leaflet provided by the

Medica Depot team for your convenience. The intellectual property rights, including but not limited to copyrights, trademarks, and other applicable rights, associated with the content of this product insert remain the sole property of the original manufacturer.

Medica Depot disclaims any responsibility for any inaccuracies, omissions, or discrepancies in this document. Your use of this product insert is subject to our terms and conditions, available on our website. By accessing and using this document, you acknowledge and agree to the mentioned terms and conditions.

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BioScience GmbH